


Product Datasheet

GFI1 antibody (orb1316176)

Description	Mouse monoclonal antibody to GFI1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GFI1 (NP_005254) produced in E.coli.
Target	GFI1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	45.1 kDa
Uniprot ID	Q99684
NCBI	NM_005263
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...