


Product Datasheet

Cytokeratin 8 (KRT8) antibody (orb1316246)

Description	Mouse monoclonal antibody to Cytokeratin 8 (KRT8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 91-381 of human KRT8(NP_002264) produced in E.coli.
Target	KRT8
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	53.5 kDa
Uniprot ID	P05787
NCBI	NM_002273
Dilution Range	IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...