


Product Datasheet

CITED1 antibody (orb1316247)

Description	Mouse monoclonal antibody to CITED1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CITED1 (NP_004134) produced in E.coli.
Target	CITED1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	19.7 kDa
Uniprot ID	Q99966
NCBI	NM_004143
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...