


Product Datasheet

MMP3 antibody (orb1316249)

Description	Mouse monoclonal antibody to MMP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 246-477 of human MMP3(NP_002413) produced in E.coli.
Target	MMP3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.2 kDa
Uniprot ID	P08254
NCBI	NM_002422
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...