


## Product Datasheet

### c Fos (FOS) antibody (orb1316250)

|                            |                                                                                    |
|----------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to c Fos (FOS)                                           |
| <b>Species/Host</b>        | Mouse                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                       |
| <b>Tested Applications</b> | IHC, WB                                                                            |
| <b>Immunogen</b>           | Full length human recombinant protein of human FOS (NP_005243) produced in E.coli. |
| <b>Target</b>              | FOS                                                                                |
| <b>Note</b>                | For research use only                                                              |
| <b>Isotype</b>             | IgG1                                                                               |
| <b>Clonality</b>           | Monoclonal                                                                         |
| <b>MW</b>                  | 40.5 kDa                                                                           |
| <b>Uniprot ID</b>          | <b>P01100</b>                                                                      |
| <b>NCBI</b>                | <b>NM_005252</b>                                                                   |
| <b>Dilution Range</b>      | IHC 1:150                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                    |



Western blot analysis of HEK293T cells w...