


Product Datasheet

ABAT antibody (orb1316252)

Description	Mouse monoclonal antibody to ABAT
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 29-323 of human ABAT(NP_065737) produced in E.coli.
Target	ABAT
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	53.2 kDa
Uniprot ID	P80404
NCBI	NM_020686
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...