


Product Datasheet

ROR beta (RORB) antibody (orb1316263)

Description	Mouse monoclonal antibody to ROR beta (RORB)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-260 of human RORB(NP_008845) produced in E.coli.
Target	RORB
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	51.9 kDa
Uniprot ID	Q92753
NCBI	NM_006914
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...