


Product Datasheet

c Fos (FOS) antibody (orb1316313)

Description	Mouse monoclonal antibody to c Fos (FOS)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Full length human recombinant protein of human FOS (NP_005243) produced in E.coli.
Target	FOS
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	40.5 kDa
Uniprot ID	P01100
NCBI	NM_005252
Dilution Range	IHC 1:250
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...