


Product Datasheet

BCL2 antibody (orb1316406)

Description	Mouse monoclonal antibody to BCL2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human BCL2 (NP_000624) produced in HEK293T cell.
Target	BCL2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P10415
NCBI	NM_000633
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...