


Product Datasheet

Desmuslin (SYNM) antibody (orb1316431)

Description	Mouse monoclonal antibody to Desmuslin (SYNM)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 329-580 of human SYNM(NP_663780) produced in E.coli.
Target	SYNM
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	172.7 kDa
Uniprot ID	O15061
NCBI	NM_145728
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...