


Product Datasheet

WIBG (PYM1) antibody (orb1316447)

Description	Mouse monoclonal antibody to WIBG (PYM1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human WIBG (NP_115721) produced in HEK293T cell.
Target	PYM1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22.5 kDa
Uniprot ID	Q9BRP8
NCBI	NM_032345
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...