

Product Datasheet

TNFRSF18 antibody (orb1316490)

Description	Mouse monoclonal antibody to TNFRSF18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 26-162 of human TNFRSF18(NP_004186) produced in E.coli.
Target	TNFRSF18
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23.5 kDa
Uniprot ID	Q9Y5U5
NCBI	NM_004195
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

