


Product Datasheet

CSNK1G1 antibody (orb1316508)

Description	Mouse monoclonal antibody to CSNK1G1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 163-422 of human CSNK1G1(NP_071331) produced in E.coli.
Target	CSNK1G1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q9HCP0
NCBI	NM_022048
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...