


Product Datasheet

PCBP1 antibody (orb1316520)

Description	Mouse monoclonal antibody to PCBP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-225 of human PCBP1(NP_006187) produced in E.coli.
Target	PCBP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37.3 kDa
Uniprot ID	Q15365
NCBI	NM_006196
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...