


Product Datasheet

HOMER2 antibody (orb1316532)

Description	Mouse monoclonal antibody to HOMER2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human HOMER2 (NP_004830) produced in HEK293T cell.
Target	HOMER2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	39.3 kDa
Uniprot ID	Q9NSB8
NCBI	NM_004839
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...