

Product Datasheet

Glucocorticoid Receptor (NR3C1) antibody (orb1316579)

Description	Mouse monoclonal antibody to Glucocorticoid Receptor (NR3C1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-265 of human NR3C1(NP_000167) produced in E.coli.
Target	NR3C1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	85.5 kDa
Uniprot ID	P04150
NCBI	NM_000176
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

