


Product Datasheet

MRRF antibody (orb1316627)

Description	Mouse monoclonal antibody to MRRF
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 56-262 of human MRRF(NP_620132) produced in E.coli.
Target	MRRF
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.1 kDa
Uniprot ID	Q96E11
NCBI	NM_138777
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...