


Product Datasheet

GRK4 antibody (orb1316633)

Description	Mouse monoclonal antibody to GRK4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-245 of human GRK4(NP_892027) produced in E.coli.
Target	GRK4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	66.4 kDa
Uniprot ID	P32298
NCBI	NM_182982
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...