


Product Datasheet

Factor I (CFI) antibody (orb1316656)

Description	Mouse monoclonal antibody to Factor I (CFI)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 340-583 of human CFI(NP_000195) produced in E.coli.
Target	CFI
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	63.4 kDa
Uniprot ID	P05156
NCBI	NM_000204
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...