


Product Datasheet

CLINT1 antibody (orb1316668)

Description	Mouse monoclonal antibody to CLINT1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-290 of human CLINT1(NP_055481) produced in E.coli.
Target	CLINT1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q14677
NCBI	NM_014666
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...