


Product Datasheet

SUOX antibody (orb1316687)

Description	Mouse monoclonal antibody to SUOX
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 320-545 of human SUOX(NP_000447) produced in E.coli.
Target	SUOX
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	60.1 kDa
Uniprot ID	P51687
NCBI	NM_000456
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...