


Product Datasheet

PROCR antibody (orb1316756)

Description	Mouse monoclonal antibody to PROCR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-210 of human PROCR (NP_006395) produced in HEK293T cell.
Target	PROCR
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	24.7 kDa
Uniprot ID	Q9UNN8
NCBI	NM_006404
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...