


Product Datasheet

PAGE1 antibody (orb1316801)

Description	Mouse monoclonal antibody to PAGE1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PAGE1 (NP_003776) produced in E.coli.
Target	PAGE1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	16 kDa
Uniprot ID	O75459
NCBI	NM_003785
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...