


Product Datasheet

hnRNP L (HNRNPL) antibody (orb1316828)

Description	Mouse monoclonal antibody to hnRNP L (HNRNPL)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 270-502 of human HNRNPL(NP_001524) produced in E.coli.
Target	HNRNPL
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64 kDa
Uniprot ID	P14866
NCBI	NM_001533
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...