

Product Datasheet

PIM2 antibody (orb131688)

Description

Goat polyclonal to PIM2

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

PIM2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

In transfected HEK293 transiently expressing Human PIM2 (myc and DYKDDDDK tagged), a band of approx. 38kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration, 0.5-1 µg/ml.WB: In transfected HEK293 transiently expressing Human PIM2 (myc and DYKDDDDK tagged), a band of approx. 38kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration, 0.5-1 µg/ml.

Protein Sequence

QTPAEDVPLNPSK

Clonality

Polyclonal

MW

34.2

Entrez


11040

Dilution Range


ELISA: 1:16000, WB: 0.5-1 µg/ml, IHC-P: 5 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of HEK293 lysate (...)



orb131688 Flow cytometric analysis of pa...



orb131688 Immunofluorescence analysis of...