

Product Datasheet

Estrogen Related Receptor beta (ESRRB) antibody (orb1316880)

Description Mouse monoclonal antibody to Estrogen Related Receptor beta (ESRRB)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 1-252 of human ESRRB (NP_004443) produced in E.coli.

Target ESRRB

Note For research use only

Isotype IgG2b

Clonality Monoclonal

Uniprot ID [O95718](#)

NCBI [NM_004452](#)

Dilution Range WB 1:500

Expiration Date 12 months from date of receipt.

