


Product Datasheet

PLAT antibody (orb1316884)

Description	Mouse monoclonal antibody to PLAT
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 157-426 of human PLAT(NP_000921) produced in E.coli.
Target	PLAT
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	60.5 kDa
Uniprot ID	P00750
NCBI	NM_000930
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...