


Product Datasheet

p53 (TP53) antibody (orb1316897)

Description	Mouse monoclonal antibody to p53 (TP53)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TP53 (NP_000537) produced in SF9 cell.
Target	TP53
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43.5 kDa
Uniprot ID	P04637
NCBI	NM_000546
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...