

Product Datasheet

ASPA antibody (orb1316934)

Description	Mouse monoclonal antibody to ASPA
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 77-313 of human ASPA (NP_000040) produced in E.coli.
Target	ASPA
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	35.6 kDa
Uniprot ID	P45381
NCBI	NM_000049
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...