


## Product Datasheet

# Constitutive androstane receptor (NR1I3) antibody (orb1316943)

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Constitutive androstane receptor (NR1I3)                                                   |
| <b>Species/Host</b>        | Mouse                                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                            |
| <b>Tested Applications</b> | IHC                                                                                                                     |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 1-236 of human NR1I3 (NP_001070950) produced in E.coli. |
| <b>Target</b>              | NR1I3                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                   |
| <b>Isotype</b>             | IgG2a                                                                                                                   |
| <b>Clonality</b>           | Monoclonal                                                                                                              |
| <b>MW</b>                  | 40.3 kDa                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q14994</a>                                                                                                  |
| <b>NCBI</b>                | <a href="#">NM_001077482</a>                                                                                            |
| <b>Dilution Range</b>      | IHC 1:150                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                         |



Immunohistochemical staining of paraffin...