


## Product Datasheet

### Mesothelin (MSLN) antibody (orb1316973)

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Mesothelin (MSLN)                                                                        |
| <b>Species/Host</b>        | Mouse                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Tested Applications</b> | WB, xMAP                                                                                                              |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 270-493 of human MSLN (NP_005814) produced in E.coli. |
| <b>Target</b>              | MSLN                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                 |
| <b>Isotype</b>             | IgG1                                                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                                            |
| <b>MW</b>                  | 62.7 kDa                                                                                                              |
| <b>Uniprot ID</b>          | <a href="#">Q13421</a>                                                                                                |
| <b>NCBI</b>                | <a href="#">NM_005823</a>                                                                                             |
| <b>Dilution Range</b>      | WB 1:2000                                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                       |



Western blot analysis of HEK293T cells w...