


Product Datasheet

SAE1 antibody (orb1317020)

Description	Mouse monoclonal antibody to SAE1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SAE1 (NP_005491) produced in HEK293T cell.
Target	SAE1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	38.3 kDa
Uniprot ID	Q9UBE0
NCBI	NM_005500
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...