


Product Datasheet

HIP55 (DBNL) antibody (orb1317051)

Description	Mouse monoclonal antibody to HIP55 (DBNL)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Full length human recombinant protein of human DBNL (NP_054782) produced in HEK293T cell.
Target	DBNL
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	48.1 kDa
Uniprot ID	Q9UJU6
NCBI	NM_014063
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...