


Product Datasheet

CIZ1 antibody (orb131706)

Description	Goat polyclonal to CIZ1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CIZ1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 90kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 91.2kDa according to NP_001124490.1). Recommended concentration: 0.3-1 µg/ml.WB: Approx 90kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 91.2kDa according to NP_001124490.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QPNTQDKTPSKVTAR
Clonality	Polyclonal
MW	100; 93.6; 93.1; 91.2; 106; 88.5
Entrez	25792
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa nuclear lysate...