


Product Datasheet

DMTF1 antibody (orb131708)

Description	Goat polyclonal to DMTF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DMTF1 (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 70kDa band observed in Human Lung and Testis lysates (calculated MW of 74.9kDa according to NP_001135798.1). Recommended concentration: 1-3 µg/ml.WB: Approx 70kDa band observed in Human Lung and Testis lysates (calculated MW of 74.9kDa according to NP_001135798.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	CPQNEADEIDSEDS
Clonality	Polyclonal
MW	84.5; 74.9
Entrez	9988
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung lysa...