


Product Datasheet

MMP7 antibody (orb1317098)

Description	Mouse monoclonal antibody to MMP7
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 95-267 of human MMP7 (NP_002414) produced in E.coli.
Target	MMP7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	19.1 kDa
Uniprot ID	P09237
NCBI	NM_002423
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...