



Product Datasheet

NPHS2 antibody (orb131711)

Description	Goat polyclonal to NPHS2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NPHS2 / SRN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 45kDa band observed in Human and Mouse Kidney lysates (calculated MW of 42.2kDa according to NP_055440.1). Recommended concentration: 0.3-1 µg/ml.WB: Approx 45kDa band observed in Human and Mouse Kidney lysates (calculated MW of 42.2kDa according to NP_055440.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SPSKPVEPLNPKK
Clonality	Polyclonal
MW	42.2
Entrez	7827
Dilution Range	ELISA: 1:64000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of taining of Huma...



Western blot analysis of staining of Pig...



Western blot analysis of Human (A) and M...