


Product Datasheet

MET antibody (orb1317130)

Description	Mouse monoclonal antibody to MET
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 125-421 of human MET (NP_000236) produced in E.coli.
Target	MET
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	153 kDa
Uniprot ID	P08581
NCBI	NM_000245
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...