


Product Datasheet

hHR23b (RAD23B) antibody (orb1317142)

Description	Mouse monoclonal antibody to hHR23b (RAD23B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-253 of human RAD23B (NP_002865) produced in E.coli.
Target	RAD23B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43 kDa
Uniprot ID	P54727
NCBI	NM_002874
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...