

Product Datasheet

MALT1 antibody pairs (orb131715)

Description	Rabbit polyclonal to MALT1. MALT1 antibodies represent a pair of fully characterized antibodies that...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal or C-terminal region of human MALT1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with GST-fusion proteins containing either the N-terminal or the C-terminal region of human MALT1.
Clonality	Polyclonal
Purity	Affinity purification
MW	92
Uniprot ID	Q9UDY8
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.