


Product Datasheet

Calpain 5 (CAPN5) antibody (orb1317171)

Description	Mouse monoclonal antibody to Calpain 5 (CAPN5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 274-578 of human CAPN5 (NP_004046) produced in E.coli.
Target	CAPN5
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	73 kDa
Uniprot ID	O15484
NCBI	NM_004055
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...