


Product Datasheet

XRCC3 antibody (orb1317190)

Description	Mouse monoclonal antibody to XRCC3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-235 of human XRCC3 (NP_005423) produced in E.coli.
Target	XRCC3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37.7 kDa
Uniprot ID	O43542
NCBI	NM_005432
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...