


Product Datasheet

Galectin 3 (LGALS3) antibody (orb1317234)

Description	Mouse monoclonal antibody to Galectin 3 (LGALS3)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human LGALS3 (NP_002297) produced in E.coli.
Target	LGALS3
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	P17931
NCBI	NM_002306
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...