


Product Datasheet

Ocular development associated gene (GATAD1) antibody (orb1317254)

Description	Mouse monoclonal antibody to Ocular development associated gene (GATAD1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 68-269 of human GATAD1 (NP_066990) produced in E.coli.
Target	GATAD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	28.5 kDa
Uniprot ID	Q8WUU5
NCBI	NM_021167
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...