

Product Datasheet

SH3RF1 antibody (orb131734)

Description	Rabbit polyclonal to SH3RF1. E3 ubiquitin-protein ligase SH3RF1 belongs to the SH3RF family. SH3RF1 ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SH3RF1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human SH3RF1.
Clonality	Polyclonal
Purity	Affinity purification
MW	93
Uniprot ID	Q7Z6J0
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.