

Product Datasheet

UHRF1 antibody pairs (orb131739)

Description	Rabbit polyclonal to UHRF1. UHRF1 antibodies represent a pair of fully characterized antibodies that...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region or the C-terminal region of human UHRF1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with GST-fusion proteins containing either the middle region or the C-terminal region of human UHRF1.
Clonality	Polyclonal
Purity	Affinity purification
MW	90
Uniprot ID	Q96T88
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.