

Product Datasheet

UHRF1 antibody (orb131741)

Description	Rabbit polyclonal to UHRF1. E3 ubiquitin-protein ligase UHRF1 is a multidomain protein that acts as ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human UHRF1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the middle region of human UHRF1.
Clonality	Polyclonal
Purity	Affinity purification
MW	90
Uniprot ID	Q96T88
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.