


Product Datasheet

TFB1M antibody (orb1317426)

Description	Mouse monoclonal antibody to TFB1M
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 116-346 of human TFB1M (NP_057104) produced in E.coli.
Target	TFB1M
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	39.4 kDa
Uniprot ID	Q8WVM0
NCBI	NM_016020
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...