

Product Datasheet

Centrin 1 (CETN1) antibody (orb1317432)

Description	Mouse monoclonal antibody to Centrin 1 (CETN1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CETN1 (NP_004057) produced in E.coli.
Target	CETN1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	19.4 kDa
Uniprot ID	Q12798
NCBI	NM_004066
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

