


Product Datasheet

C10orf63 (ENKUR) antibody (orb1317478)

Description	Mouse monoclonal antibody to C10orf63 (ENKUR)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ENKUR (NP_659447) produced in E.coli.
Target	ENKUR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.3 kDa
Uniprot ID	Q8TC29
NCBI	NM_145010
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...